



Professional Packaging Systems

2010 S. Great Southwest Parkway
Grand Prairie, Texas 75051
www.ProPac.com

Contact

Our Packaging Experts
solutions@propac.com
888-318-0083

Starview FAB-PH Medical & Pharmaceutical Blister Sealer



At a Glance

- Ideal for medium to higher volume medical and pharmaceutical packaging
- Six and eight-station designs with 14 x 18 and 18 x 24 inch sealing areas
- Process controls for validation
- Epoxy and stainless steel utilized for clean room applications
- Heavy-duty industrial design

Buy your [blisters and cards](#) from us, and get a discount on this Starview medical blister sealer!

Pricing

Call us at **888-318-0083**

Many options are available. Complete details will need to be finalized to determine requirements and final system costs. Work with your Pro Pac representative to build this machine to your specifications.

High-Volume Medical Blister Packaging

Starview FAB-PH Series rotary type heat sealers are designed for medium to higher volume production requirements where automatic feeding and unloading is desired. FAB-PH sealing machines are offered in six and eight-station designs with standard sealing areas of 14 x 18 and 18 x 24 inches. Larger sealing areas are available on request.

These medical rotary sealing machines can be configured to produce all of the medical and pharmaceutical packages including Tyvek, foils and medical kraft to PETG or HIPS trays and typical wallet packs for clinical trials or controlled dosage.

Features

A wide range of standard and custom features available to match customer requirements.

Generous use of stainless steel and powder coated steel in the structure for the most durable finish in the industry.

Color touch screen HMI for the ultimate in operator convenience and accessibility. Includes storage for 99 databases and on-screen data tracking.

A Starview exclusive remote access via standard Ethernet hardwire or optional Wi-Fi router gives management access to real-time process information from the machine.

On-screen temperature calibration standard with optional on-screen pressure calibration available to simplify validation using direct access validation ports.

Precision indexing rotary table for precise sealing fixture to matching relieved face plate alignment.

Ultra-high pressure sealing presses used in conjunction with matched aluminum tooling to give consistent sealing patterns with any heat sealable medical lidding materials. Hydro-pneumatic press available for cold seal materials.

Industry leading design, construction, features and customer support.

Please [contact Pro Pac](#) or call 888-318-0083 for your blister sealing equipment.

Specifications of the FAB-PH Models

Specifications	FAB/PH6-1418	FAB/PH6-1824	FAB/PH8-1418	FAB/PH8-1824
Sealing Area	14" x 18"	18" x 24"	14" x 18"	18" x 24"
Open Stations	2	2	4	4
Table Index	Automatic	Automatic	Automatic	Automatic
Heating System	2 Zones	2 Zones	2 Zones	2 Zones
Automatic Feeders	2	2	2	2
Finished Package Unloader	Standard	Standard	Standard	Standard
Maximum Blister Tray Depth	3"	3"	3"	3"
Cycle Speed	1 - 18 CPM	1 - 14 CPM	1 - 18 CPM	1 - 16 CPM
Electrical	208-230V, 3Ph, 25Amp	208-230V, 3Ph, 45Amp	208-230V, 3Ph, 25Amp	208-230V, 3Ph, 45Amp