



## Professional Packaging Systems

2010 S. Great Southwest Parkway  
Grand Prairie, Texas 75051  
[www.ProPac.com](http://www.ProPac.com)

## Contact

Our Packaging Experts  
[solutions@propac.com](mailto:solutions@propac.com)  
888-318-0083

# Starview MR/PH2 Medical & Pharmaceutical Blister Sealer



## At a Glance

- Ideal for medium-volume medical and pharmaceutical packaging
- High-pressure sealing presses for consistent results with Tyvek lidstock
- Process controls for validation
- Epoxy and stainless steel utilized for clean room applications
- Heavy-duty industrial design

Buy your [blisters and cards](#) from us, and get a discount on this Starview medical blister sealer!

## Pricing

Call us at **888-318-0083**

Many options are available. Complete details will need to be finalized to determine requirements and final system costs. Work with your Pro Pac representative to build this machine to your specifications.

## Quickly Blister Pack Medical Products in a Small Foot Print

Starview MR/PH2 Series two-station manual rotary type heat sealers are designed for medium-volume production requirements. Machines have standard sealing areas of 14 x 18 and 18 x 24 inches. These manual rotary sealing machines can be configured to produce all of the medical/pharmaceutical packages including Tyvek, foils and medical kraft to PETG or HIPS trays and typical wallet packs for clinical trials or controlled dosage.

The MR/PH2 features one sealing station with the second station used for load and unload. These high-quality sealing machines are ideal for start-up companies, back up equipment, cellular production lines, and low-volume in-house packaging.

## Features

A wide range of standard and custom features available to match customer requirements.

Generous use of stainless steel and powder coated steel in the structure for the most durable finish in the industry.

Color touch screen HMI for the ultimate in operator convenience and accessibility. Includes storage for 99 databases and on-screen data tracking.

A Starview exclusive remote access via standard Ethernet hardwire or optional Wi-Fi router gives management access to real-time process information from the machine.

On-screen temperature calibration standard with optional on-screen pressure calibration available to simplify validation using direct access validation ports.

Compact manual rotary table design conserves clean room space.  
 Ultra-high pressure sealing presses used in conjunction with matched aluminum tooling to give consistent sealing patterns with any heat sealable medical lidding materials.  
 Industry leading design, construction, features and customer support.

Please [contact Pro Pac](#) or call 888-318-0083 for your blister sealing equipment.

## Specifications of the MR/PH2 Models

Specifications	MR/PH2-1418	MR/PH2-1824
Sealing Area	14" x 18"	18" x 24"
Open Stations	1	1
Heating System	2 Zones	2 Zones
Table Index	Manual	Manual
Maximum Blister Tray Depth	3"	6"
Cycle Speed	1 - 6 CPM	1 - 6 CPM
Electrical	208-230V, 1Ph, 45Amp	208-230V, 1Ph, 25Amp or 208-230V, 3Ph, 25Amp